A-1341 Room G, 10/17/2000 2:00 PM - 4:00 PM (PS) Interaction between Volatile Anesthetics and Hypoxia in Porcine Tracheal Smooth Muscle Michiaki Yamakage, MD, PbD; Xiangdong Chen, MD; Naoki Tsujiguchi, MD; Yasuhiro Kamada, MD; Akiyoshi Namiki, MD, PbD, Anesthesiology, Sapporo Medical University School of Medicine, Sapporo, Hokkaido, Japan. Hypoxia inhibits airway smooth muscle contraction independently of intracellular Ca²⁺.

A-1342 Room G, 10/17/2000 2:00 PM - 4:00 PM (PS) Both PKA and PKG are Essential in Regulating Basal Ciliary Beat Frequency in Rat Tracheal Epithelial Cells Xinhua Zhan, M.D., Ph.D; Dechun Li, M.D., Ph.D; Roger A. Johns, M.D., Department of Anesthesiology, Union Affiliated Hospital of Tongji Medical University, Wuhan, Hubei, China. Inhibition of PKA or PKG did not change CBF. Inhibition of both decreased CBF significantly.

Respiration: Airway / Pulmonary Vascular Smooth Muscle

A-1343 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Mechanisms of Bronchoprotection by Propofol In Vivo: Effects of Preservative Robert H. Brown, M.D., M.P.H.; Robert S. Greenberg, M.D.; Elizabeth M. Wagner, Ph.D., Anesthesiology and Critical Care Medicine, Johns Hopkins Medical Institutions, Baltimore, M.D., Metabisufite (MBS) increases direct and neurally-mediated bronchoconstriction. Propofol without (MBS) decreased neurally-mediated bronchoconstriction.

A-1344 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Intraveous Colforsin Daropate Prevents Thiamylal-Fentanyl-Induced Bronchoconstriction Zen'ichiro Wajima, MD,PhD; Tatsusuke Yoshikawa, MD,PhD; Akira Ogura, MD,PhD; Kazuyuki Imanaga, MD; Tetsuo Inoue, MD,PhD, Dept. of Anesthesia, Chiba Hokusoh Hospital, Nippon Medical School, Inba, Chiba, Japan. After thiamylal and fentanyl, airway resistance unchanged by iv. colforsin daropate.

A-1345 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Different Inhibitory Effects of Volatile Anesthetics on T- and L-type Calcium Channels in Airway Smooth Muscle Michiaki Yamakage, MD, PbD; Xiangdong Chen, MD; Naoki Tsujiguchi, MD; Yasubiro Kamada, MD; Akiyoshi Namiki, MD, PbD, Anesthesiology, Sapporo Medical Univ. Sch. of Med., Sapporo, Hokkaido, Japan. Volatile anesthetics have more inhibitory effects on T-type calcium channels.

A-1346 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Effect of Hexanol on Calcium Sensitivity in Airway Smooth Muscle Hayashi Yoshimura, MD; Keith A. Jones, MD; William J. Perkins, MD; Shosuke Takahashi, MD; David O. Warner, MD, Department of Anesthesiology, Mayo Clinic and Foundation, Rochester, MN, United States. Hexanol inhibits agonist-induced Ca²⁺-sensitization in permeabilized airway smooth muscle

A-1347 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Propofol Increases Myofilament Calcium Sensitivity Via Protein Kinase C Activation Satoru Tanaka, MD; Izumi Kondo, MD; Derek Damron, PhD; Paul Murray, PhD, Anesthesiology Research, Cleveland Clinic Foundation, Cleveland, OH, United States. Propofol decreases [Ca²⁺]_i but increases myofilament Ca²⁺ sensitivity in pulmonary artery smooth muscle.

A-1348 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Endothelin-1 Induces Sustained Contraction without Myosin Light Chain Phosphorylation in Porcine Pulmonary Artery Chie Sakibara, M.D., Ph.D.; Tetsuya Kai, M.D., Ph.D.; Shosuke Takabashi, M.D., Ph.D., Department of Anesthesiology, Kyushu University, Fukuoka, Japan. Ca²⁺-insensitive contraction might be mediated by phosphorylation of the protein that is sensitive to H-7 and staurosporine.

A-1349 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Selective Inhibition of the Triphasic HPV Response by Inhalational Anesthetics Bryan E. Marshall, MD; Masami Ozaki, MD; James E. Baumgardner, MD PhD; Carol Marshall, PhD, Center for Anesthesia Research, University of Pennsylvania Health System, Philadelphia, PA, United States. These studies suggest that isoflurane and halothane inhibit the force sensitization mechanisms of HPV with the greatest potency.

A-1350 Room 309, 10/16/2000 2:00 PM - 3:30 PM (PD) Endogenous Nitric Oxide Does Not Play a Significant Role in the Maintenance of Basal Pulmonary Microvascular Tone Sanjay Bbatia, MB;BS; Paul Murray, PbD; David DeFily, PbD, Center for Anesthesiology Research, Cleveland Clinic Foundation, Cleveland, OH, United States. Nitric oxide does not regulate the baseline diameter of 40-140µm pulmonary arteries in vivo.

Respiration: Lung Injury / Cellular Physiology

A-1351 Room 224–226, 10/16/2000 3:30 PM - 5:00 PM (PD) Human Airway Smooth Muscle Expresses mRNA Encoding Seven Subtypes of Adenylyl Cyclase Charles W. Emala, M.D.; Dingbang Xu, Anesthesiology, Columbia University, New York, NY, United States. Messenger RNA encoding 7 of the 9 known subtypes of adenylyl cyclase, the target protein of b-adrenoceptor bronchodilating agents, was identified by RT-PCR in human airway smooth muscle.

A-1352 Room 224–226, 10/16/2000 3:30 PM - 5:00 PM (PD) Biochemical Mechanism for the Hydrogen Peroxide-Induced Inhibition of Calcium Sensitivity in Airway Smooth Muscle Robert R. Lorenz; William J. Perkins; David O. Warner; Keith A. Jones, Anesthesiology, Mayo Clinic and Foundation, Rochester, MN, United States. H₂O₂ relaxes airway smooth muscle in part by decreasing regulatory myosin light chain phosphorylaton and Ca²⁺ sensivity.

A-1353 Room 224–226, 10/16/2000 3:30 PM - 5:00 PM (PD) Nitric Oxide Production and Stimulation of Ciliary Motility in Rat Tracheal Ciliated Epithelial Cells Gotaro Shirakami, M.D.; Dechun Li, M.D., Pb.D. D; Kazubiko Fukuda, M.D.; Roger A. Johns, M.D., Department of Anesthesia, Kyoto University Hospital, Kyoto, Japan. Nitric oxide production was detected and it was positively correlated with cilialy motility in cultured rat tracheal epithelial cells.

A-1354 Room 224–226, 10/16/2000 3:30 PM - 5:00 PM (PD) Hypoxia Mediates Upregulation of ICAM-1 in Rat Alveolar Macrophages Beatrice Beck-Schimmer, M.D.; Caveb Madjdpour, M.D.; Thomas Pasch, M.D.; Ralph C. Schimmer, M.D., Institutes of Anesthesiology and Physiology, University of Zurich, Zurich, Switzerland. Hypoxia upregulates ICAM-1 expression in alveolar macrophages time-and concentration-dependent.